

# PC-66A

## Handheld Pulse Oximeter



### Features

- Measuring SpO2, Pulse Rate and Temperature(optional)
- 3.5 inch TFT display of PI%, plethysmograph and pulse bar
- Adult, pediatric and neonate measurement
- Lithium battery up to 18 hours
- 500 hours SpO2 data storage
- PC software for data management



### Technical Specifications

<b>SpO<sub>2</sub></b>	
Technique:	Dual-wavelength LED
Display range:	0%~100%
Measuring accuracy:	±3%(range from 70% to 100%)

<b>Pulse Rate</b>	
Measuring range:	30~250bpm
Measuring accuracy:	± 2bpm or ±2%

<b>Temperature(option)</b>	
Measuring range:	32℃~43℃
Measuring accuracy:	± 0.2℃
Measuring time:	≤ 5 seconds
Measuring probe:	Infrared ear probe

<b>Perfusion Index</b>	
PI display range:	0.2%~20%

<b>Measuring Mode</b>	
Continuous or Spot check	

<b>Data Storage</b>	
Up to 500 hours	

<b>Display</b>	
3.5" color TFT display with rotation	

<b>Power Supply</b>	
2000mAh lithium battery	

<b>Charging</b>	
Charging stand or USB cable	

<b>Dimension</b>	
158mm x 73mm x 25mm	

<b>NetWeight</b>	
230g	

<b>Data Upload</b>	
USB cable or wireless	

<b>Operating Environment</b>	
Operating Temperature: 5℃~40℃	
Operating Humidity: 15%-93%	



Adult clip sensor



Adult soft sensor



Pediatric clip sensor



Pediatric soft sensor



Neonate wrap sensor



Infrared ear thermometer



FDA 510(k)

Shenzhen Creative Industry Co., Ltd.  
 Floor 5, BLD 9, Baiwangxin High-Tech Industrial Park, Songbai Road, Xili Street  
 Nanshan District, 518110 Shenzhen, PEOPLE'S REPUBLIC OF CHINA  
 Tel: +86 755 2643 3514 Fax: +86 755 2643 0390  
 www.creative-sz.com info@creative-sz.com