

Oxygen Generator



M9

Oxygen Generator

Features:

- **Upgraded oxygen generation technology**
IMR-4 oxygen generation system, equipped with molecular sieve imported from France, assuring faster, thicker and purer oxygen output
- **12-inch visual oxygen regulation**
12-inch backlight screen, with a unique oxygen control system, fine-tuned to 0.1L/min
- **Voice broadcast + one-click for help**
One-button respirator for the elderly, one-button for help, making oxygen inhalation more secure
- **Long-distance wireless remote control**
Long remote control distance, 360 degrees without dead angle and of high sensitivity
- **Timed oxygen supplement for up to 500 minutes**
Timed oxygen supplement per own needs, automatic shutdown after completing oxygen supplement, which is easy to operate

Technical Parameters

Product dimension:	430x228x554mm
Net weight:	12kg
Flow range:	1-9L/min adjustable
Maximum outlet pressure:	About 47±10% kPa
Electrical classification:	Class II Type B
Running noise:	< 51dB(A)
Oxygen concentration:	93%±3%
Nominal voltage:	AC 220V 50Hz
Input power:	250VA

Ordering Information

Catalogue No.

M9